

This new and advanced RF system will help make your skin tightened, smooth and fresh, as well as improving contour of body and face.

Tip: For this machine to work properly & effectively please use [conductive gel](#) during the treatment.



How this Machine Works, in Principle:

- Technically speaking, Radio Frequency (RF) is a high frequency and alternating electromagnetic wave. This Machine uses a 5MHZ variable frequency system that rapidly alternates the electric field in the skin from positive to negative in order to set the optimum applicable frequency.
- The RF (radio frequency) produces heat (45 °C -65 °C) in the skin tissue as electrical resistance converts the electric current to thermal energy deep within the dermis. Within a short time collagen within thermally modified deep tissues shrinks, which can

tighten the flabby skin and stimulate collagen to grow. Subsequent neocollagenesis then **further tightens the dermal tissue and reduces wrinkles.** This process is what largely help to achieve the **best treatment result - skin that is tightened, smooth, and fresh.**

- The 5MHZ RF wave passes about 2 to 4cm into the fatty skin layer causing breakdown of the fat membrane, so **fat is released and picked up by the lymphatic system and disposed of** by the liver for excretion. Thus **cellulite and body contouring can be improved!**

Treatment Benefits Summary:

- Reduce wrinkles and fine lines on face and body
- Lift and contour the eye area
- remove eye bags and dark circles
- Reduce apparent pore size
- Reduce swelling of the skin
- Improve contour of the body and reduce cellulite
- Improve contour of the face, achieve V-shaped face
- Lose double chin
- Improve skin discoloration issues, reveal a rosy face color

How To Use, Basics:

* detailed instructions manual will be provided with the machine

HOW TO USE - BODY

1. Attach the larger handle to the machine
2. Plug in the power
3. Apply gel to the area
4. Adjust the mode for body and timer
5. Select level 5+ for the body
6. Massage in small circles focusing on one section
7. If the gel dries out, spray the area with water.

HOW TO USE - FACE

1. Select attachment and screw into the handle
2. Attach handle to the machine
3. Plug in the power
4. Apply gel to one side of your face
5. Adjust the mode for face and set the timer
6. Select level 3 or 4 for the face area
7. Massage in small gentle circles
8. If the gel dries out, spray the area with water.

HOW TO USE - EYES

1. Select attachment and screw into the handle
2. Attach handle to the machine
3. Plug in the power
4. Apply gel to one side of your face
5. Adjust the mode for eyes and set the timer
6. Select level 3 or 4 for the eye area
7. Massage in small gentle circles
8. If the gel dries out, spray the area with water.





Included:

- 1 x Main Machine Unit
- 1 x Black Support Stand
- 2 x Probe (one for Body, one for Face and Eyes)
- 1 x Attachment for Face
- 1 x Attachment for Eyes
- 1 x Power Cord Plug
- 1 x Instructions Manual
- 1 x Black Storage Box









**4 cores connector
diameter: 0.9cm**



with red light



4.8cm/1.88"

Body RF Wand

G.W.: 0.25kg



Face RF tip

G.W: 0.055kg

(no red light)



1.6cm/0.62"

Eye area RF tip

G.W: 0.04kg

(no red light)